

# Warning:

Wearing gloves is highly recommended when handling the kit contents.

*GeBAflex-tube* is covered by the WO0190731 patent application assigned to Gene Bio-Application Ltd.

*GeBAflex-tubes* are autoclaved and are bacterial free.

All kit buffers are filtered, autoclaved and are bacterial free.

***GeBAflex-tubes* membrane is ultra-clean, sulfur and heavy metal free and EDTA treated.**

# ***Mini GeBAflex-tube (250 µl) Gel Extraction and Dialysis Kit Handbook***

## **Table of Contents**

<b>Applications</b> .....	4
<b>Mini <i>GeBAflex-tube</i> kit Contents</b> .....	4
<b>Storage Conditions</b> .....	5
<b>Product use Limitations</b> .....	5
<b>Quality Control</b> .....	5
<b><i>GeBAflex-tube</i></b> .....	5
<b>Yield of Molecule Recovery</b> .....	6
<b>Specifications</b> .....	6
<b>Protein Extraction from Polyacrylamide Gel with <i>GeBAflex-tube</i></b> ..	7
Procedure.....	7
Elution time table .....	10
Protein precipitation protocols.....	11
Trichloroacetic acid (TCA) precipitation procedure .....	11
MS precipitation procedure.....	12
<b>Protein Extraction from SDS-PAGE compatible with MALDI-MS with <i>GeBAflex-tube</i></b> .....	13
Introduction .....	13
Procedure.....	14
Protein Precipitation protocol for analysis with MALDI-MS.....	16
<b>Dialysis with the <i>GeBAflex-tube</i></b> .....	17
Procedure.....	17

<b>Sample Concentration by Evaporation with <i>GeBAflex-tube</i>...</b>	.19
<b>DNA and RNA Extraction from Gel with <i>GeBAflex-tube</i> .....</b>	19
Procedure.....	20
Elution timetables .....	21
DNA or RNA precipitation .....	23
<b>Troubleshooting Guide .....</b>	24
<b>Ordering Information .....</b>	25

## Applications

- Extraction of proteins, RNA, DNA or oligonucleotides (>20 nt) from polyacrylamide, agrose or any gel matrix in any running buffer.
- Extraction of protein-protein, DNA-protein or RNA-protein complexes.
- Dialysis or buffer exchange of volumes between 10-250  $\mu$ l
- Preparation of protein samples for MALDI-MS.
- 2D gel dots extraction.
- Samples concentration

## Mini *GeBAflex-tube* Kit Contents

Mini <i>GeBAflex-tube</i> kit (electro-elution & dialysis)	For 10	For 30
Mini <i>GeBAflex-tubes</i> 250 $\mu$ l	10 units	30 units
Supporting tray (for electro elution protocols)	1	1
Floating rack (for dialysis protocol)	1	1
MS Buffer	250 $\mu$ l	750 $\mu$ l
20 % TCA Buffer	2.5 ml	7.5 ml
3M potassium acetate, pH 5.2 (KAc) Buffer	300 $\mu$ l	900 $\mu$ l
Information and protocols handbook	1	1

\* **Note:** Kits for **dialysis do not** contain Supporting tray and Buffers.

## Storage Conditions

*GeBAflex-tube* kit must be stored in a dry place at room temperature (15-25°C). Under these conditions, *GeBAflex-tube* kit can be stored for up to 12 months without any deterioration in performance and quality. For longer storage time, it is recommended that the *GeBAflex-tube* kit be stored in a cool place (refrigerator), at relative humidity of 35% at least.

## Product Use Limitations

*GeBAflex-tube* kit is developed, designed and sold for research purposes only. It is not to be used for human diagnostic purposes or drug production nor for producing any substance intended to be administered to humans unless expressly cleared for that purpose by the Food and Drug Administration in the USA or the appropriate regulatory authorities in the country of use. All due care and attention should be exercised in the handling of materials described in this text.

## Quality Control

The performance of *GeBAflex-tube* kit is regularly monitored. *GeBAflex-tube* kit is tested by using it for extraction of Proteins, DNA and RNA fragments of various sizes from either agarose or polyacrylamide gel. *GeBAflex-tube* kit is tested also for simple dialysis of salts or buffer exchange. The quality of the isolated Protein, DNA and RNA fragments or of the sample after dialysis is checked by several assays commonly used for proteins, nucleic acids and dialysis. Determining the recovery from a specific amount of loaded sample tests the quality and efficiency of the *GeBAflex-tube* membrane.

## *GeBAflex-tube*

The device combines two modes of action, electro-elution of macromolecules from polyacrylamide or agarose gel and dialysis or buffer exchange at volume samples between 10-250 µl. This device allows rapid and high performance at either mode and extracts the macromolecules without any contamination.

## Yield of Molecule Recovery

DNA or RNA from agarose gel	90%
DNA or RNA from polyacrylamide gel	90%

Protein from SDS-PAGE	70%
-----------------------	-----

## Specifications

Regenerated cellulose membrane cut-off	6000- 8000 or 12000 14000 MWCO
Membrane	Ultra-clean Sulfur and heavy metal free EDTA treated
Tube maximum volume capacity	250 $\mu$ l
Minimum amount of protein at the start of extraction	0.5 $\mu$ g
Maximum size of the gel slice that can be inserted into the tube	0.4 cm x 1.1 cm
Volume of sample for dialysis	10-250 $\mu$ l

## **Protein Extraction from Polyacrylamide Gel with *GeBAflex-tube***

**IMPORTANT:** Fixation of proteins before electro elution (**e.g. fixation with methanol, acetic acid, etc**) is not recommended; fixation greatly reduces extraction yield. A sensitive protein staining solution, **SeeBand (from Gene Bio-Application Ltd, see Ordering Information, page 30)**, is a good staining reagents, as it permanently stain the gel without undue fixing of the protein.

## Procedure

1. Fill the *GeBAflex-tube* with 250  $\mu\text{l}$  of  $\text{dH}_2\text{O}$ ; incubate for at least 5 min. empty the tube.

**IMPORTANT:** Check carefully that no  $\text{dH}_2\text{O}$  is leaking from the tube. Absorbent of water, by the dry membrane, cause the decrease in water level.

2. After staining the gel, excise the gel slice containing the protein with a clean, sharp scalpel.

Minimize the size of the gel slice by removing extra gel. Maximum gel slice size 0.4 cm x 1.1 cm.

Using the *SeeBand* staining solution (see Ordering Information, page 30) will result in highest recovery yield of proteins from the gel.

3. Transfer the gel slice to a *GeBAflex-tube*. Fill the tube with protein-running buffer (about 250  $\mu\text{l}$ ). Close the tube gently with the cup (do not apply force).

**IMPORTANT:** Avoid air bubbles in the tube. **Do not fill the tube with several gel slices**, for larger gel slices use more than one tube.

4. Place the *GeBAflex-tube* in the provided supporting tray (see Figure 1).

The supporting tray can hold 1-4 *GeBAflex-tube(s)*.


**IMPORTANT:** The two membranes of the *GeBAflex-tube* must be in perpendicular to the electric field to permit the electric current to pass through the tube.



**Figure 1:** Insertion of the *GeBAflex-tube* in the provided supporting tray.

- 5.** Place the supporting tray containing the *GeBAflex-tube(s)* in a horizontal electrophoresis tank containing protein running buffer (see Figure 2).

**IMPORTANT:** Immerse fully the *GeBAflex-tube(s)* with the tray in the buffer.



**Figure 2:** Supporting tray-containing 4 Mini *GeBAflex-tubes* in a horizontal electrophoresis tank. The arrow on the cap is positioned face up and the two membranes of the Mini *GeBAflex-tube* are in perpendicular to the electric field.

**6. Pass electric current (usually at 100 volt) until the protein exits from the gel slice.**

Electro-elution time is to be adjusted for each individual sample. It takes at least 70-80 min for BSA protein to be electro-eluted from a 10% SDS-PAGE slice, size 0.4 cm x 1.1 cm (see Table 1, page 11).

**7. Reverse the polarity of the electric current for 120 seconds.**

This step will release the protein from the membrane.

**8. Open the *GeBAflex-tube* gently, pipetting the protein-containing solution up and down carefully (at least 5 times) and transfer the solution to a clean tube.**

Do the pipetting on the inner side of the membrane.

**Important Notes:**

- i. Use the extracted protein directly.
- ii. Concentrate the extracted protein by *ProteoConN* or *ProteoConD* kits (see Ordering Information Cat # PDT030, PDT035, PNT030 and PNT035 page 28-30).
- iii. Precipitate the extracted protein by standard precipitation protocols (see page 11-12).
- iv. Dialyze directly the extracted protein with a clean *GeBAflex-tube* (see page 17).

**Elution Time Table**

The elution time depends on the size of the protein molecule to be eluted, the applied voltage, the size of gel slice, the ratio of the polyacrylamide:bisacrylamide, the intensity of the fixation of the protein to the gel by the staining solution and the

percentage of the polyacrylamide gel. **Electro-elution time at the elution step was to be adjusted for each individual sample.**

**Table 1:** Minimum time needed to extract different-sized proteins from 10% SDS-polyacrylamide gel (37.5:1 polyacrylamide:bisacrylamide) at 100 V, in 1XPRB: 0.192M Glycine, 0.025M Tris-base and 0.1% SDS.

**IMPORTANT:** proteins were visualized by staining with SeeBand cat # SB010 (see Ordering Information, page 30).

Protein (kDa)	Time (min)	Protein (kDa)	Time (min)
14.4	30-40	66.2	70-80
18.4	35-45	116	90-100
25	40-50		
35	50-60		
45	55-65		

## Protein precipitation protocols

### Trichloroacetic acid (TCA) precipitation procedure

1. Add equal volume of 20% TCA to the tube containing the extracted protein solution and mix properly.

For example, add 250  $\mu$ l of 20% TCA to a 250  $\mu$ l sample.

2. Incubate 60 min in 4°C.
3. Spin the tube at 4°C for 30 min at 14,000 RPM.
4. Discard supernatant carefully.
5. Add 0.5 ml cold acetone.

6. Incubate at  $-20^{\circ}\text{C}$  for 60 min and centrifuge the sample at  $4^{\circ}\text{C}$  for 30 min at 14,000 RPM.

To increase protein precipitation yield incubate the samples over night at  $-20^{\circ}\text{C}$ .

7. Discard supernatant and air-dry the pellet.
8. Resuspend the pellet using 0.1M NaOH (**use at least 0.1 ml to perform resuspension**).

### **MS precipitation procedure (recommended when protein-bound SDS need to be removed)**

1. Add 1:10 by volume of MS buffer to the protein containing solution and mix properly.

For example, add 25  $\mu\text{l}$  of MS buffer to a 250  $\mu\text{l}$  sample.

2. Incubate for 15 min at room temperature.
3. Add 1: 2 by volume of 20% TCA and mix properly.

For example, add 138  $\mu\text{l}$  of 20% TCA to a 275  $\mu\text{l}$  sample.

4. Incubate for 1 hour at  $4^{\circ}\text{C}$ .
5. Centrifuge the sample at  $4^{\circ}\text{C}$  for 30 min at 14,000 RPM.
6. Carefully descent the supernatant without disturbing the pellet.
7. Add 0.5 ml of ice-cold acetone.
8. Incubate at  $-20^{\circ}\text{C}$  for 30 min and centrifuge the sample at  $4^{\circ}\text{C}$  for 30 min at 14,000 RPM.

To increase protein precipitation yield incubate the samples over night at  $-20^{\circ}\text{C}$ .

9. Carefully descent the supernatant without disturbing the pellet. Air-dry the pellet.
10. Resuspend the pellet in a suitable buffer solution or 0.1M NaOH (**use at least 20  $\mu\text{l}$  to perform resuspension**).

# **Protein Extraction from Polyacrylamide Gel compatible with Matrix-assisted Laser Desorption/Ionization Mass Spectrometry (MALDI-MS) by Mini *GeBAflex-tube***

## **Introduction**

Sodium dodecylsulfate-polyacrylamide gel electrophoresis (SDS-PAGE) is a widely used technique for the separation and molecular weight estimation of individual proteins. However, the accuracy of this molecular weight determination is often inadequate for protein characterization. More recently Matrix-assisted Laser Desorption/Ionization Mass Spectrometry (MALDI-TOF-MS) has found widespread use for the determination of molecular mass of intact proteins isolated from gels. The isolation of proteins from gels with the newly developed *GeBAflex-tube* electro-elution system provides 70% recovery yields. This combination of SDS-PAGE, *GeBAflex-tube* electro-elution system and MALDI-TOF-MS is attractive. It provides a much more accurate determination of protein molecular weight. Moreover, even difficult proteins to analyze such as integral membrane proteins (hydrophobic) or high molecular mass proteins can be analyzed. This unique method provides a powerful means for characterizing endogenous proteins of wide molecular weight range separated by SDS-PAGE.

The combination of the three methods provides significantly improved protein yield and SDS free samples. The end result is a MALDI-MS analysis with greater sensitivity. The *GeBAflex-tube* tool provides high protein yield recovery, and the MS buffer contained in the *GeBAflex-tube* kit thoroughly removes the SDS.

## **Procedure**

- 1. Fill the *GeBAflex-tube* with 250  $\mu$ l of dH<sub>2</sub>O; incubate for at least 5 min., empty the tube.**

**IMPORTANT:** Check carefully that no dH<sub>2</sub>O is leaking from the tube. Absorbent of water, by the dry membrane, cause to the decrease in water level.

- 2. After staining the gel, excise the gel slice containing the protein with a clean, sharp scalpel.**

Minimize the size of the gel slice by removing extra gel.  
Maximum gel slice size 0.4 cm x 1.1 cm.

Using *SeeBand* staining solution will result in highest recovery yield of proteins from gel.

3. **Transfer the gel slice to a *GeBAflex-tube*. Fill the tube with protein running buffer: 250 mM Tricine pH 8.5, 0.025% SDS and 25 mM Tris-Base. Close the tube gently with Mini *GeBAflex-tube* cap (do not apply force).**

**IMPORTANT:** Avoid air bubbles in the tube. **Do not fill the tube with several gel slices**, for large gel slices use more than one tube.

4. **Place the *GeBAflex-tube* in the provided supporting tray (see Figure 1, page 8).**

The supporting tray can hold 1-4 *GeBAflex-tube(s)*.

5. **Place the supporting tray containing the *GeBAflex-tube(s)* in a horizontal electrophoresis tank filled with protein-running buffer: 250 mM Tricine pH 8.5, 0.025% SDS and 25 mM Tris-Base (see Figure 2 page 9).**

**IMPORTANT:** Immerse **fully** the *GeBAflex-tube(s)* with the tray in the buffer.

6. **Pass electric current at 150 volt until the protein exits from the gel slice.**

**Electro-elution time at the elution step was to be adjusted for each individual sample.** It takes at least 90-100 min for BSA protein to be electro-eluted from a 10% SDS-PAGE gel slice in the size of 0.4 cm x 1.1 cm.

For other proteins from BSA, increase electro elution time presented in **Table 1** page 11, by 30%.

The elution time depends on the size of the protein molecule to be eluted, the applied voltage, the size of gel slice, the ratio of the polyacrylamide:bisacrylamide, the intensity of the fixation of the protein to the gel by the staining solution and the percentage of the polyacrylamide gel.

7. **Reverse the polarity of the electric current for 120 seconds.**

This step will release the protein from the membrane.

8. Open the *GeBAflex-tube* gently, pipetting the protein-containing solution up and down carefully (at least 5 times) and transfer the solution to a clean tube.

Do the pipetting on the inner side surface of the membrane.

**Important notes:** See page 10.

### Protein precipitation protocol for analysis by MALDI-MS

1. Add 1:10 by volume of MS buffer to the protein containing solution and mix properly.

For example, add 25  $\mu$ l of MS buffer to a 250  $\mu$ l sample.

2. Incubate for 15 min at room temperature.
3. Add 1: 5 by volume of **50% TCA** (not supplied in the kit) and mix properly.

For example, add 55  $\mu$ l of **50% TCA** to a 275  $\mu$ l sample.

4. Incubate for 1 hour at 4°C.
5. Centrifuge the sample at 4°C for 30 min at 14,000 RPM.
6. Carefully descent the supernatant without disturbing the pellet.
7. Add 0.5 ml of ice-cold acetone.
8. Incubate at -20°C for 30 min and centrifuge the sample at 4°C for 30 min at 14,000 RPM.

To increase protein precipitation yield incubate the samples over night at -20°C.

9. Carefully descent the supernatant without disturbing the pellet. Air-dry the pellet.
10. For mass spectrometric analysis resuspend the pellet in appropriate solution compatible with MALDI-MS (protein characteristic is important for determination the appropriate solution) followed by essential dilution step according to the protocols compatible with MALDI-MS. Use at least 20  $\mu$ l to perform resuspension.

## Dialysis with Mini *GeBAflex-tube*



**Figure 3:** Dialysis with Mini *GeBAflex-tube*.

### **Procedure**

- 1. Fill the *GeBAflex-tube* with 250  $\mu\text{l}$  of  $\text{dH}_2\text{O}$ ; incubate for at least 5 min. empty the tube.**

**IMPORTANT:** Check carefully that there is no  $\text{dH}_2\text{O}$  leaking from the tube. Absorbent of water, by the dry membrane, cause to the decrease in water level.

- 2. Load sample into the *GeBAflex-tube*. Close the tube with the provided caps (do not apply force).**

Sample volume should be in the range of 10-250  $\mu\text{l}$ . If small volume is used (e.g., 10  $\mu\text{l}$ ), load the sample close to the inner membrane.

- 3. Place the loaded *GeBAflex-tube* in the supplied floating rack in a stirred beaker containing large volume (usually 100 to 1000-fold that of the sample) of the desired buffer.**

The floating rack can hold 1-4 *GeBAflex-tube(s)*.

Adjust the stir bar speed. Allow at least 30 min for each 0.1 ml of sample. Low-molecular weight salts and buffers (e.g., Tris•Cl and KPO<sub>4</sub>) equilibrate within 3 hours. Equilibration times for viscous samples will be longer.

**IMPORTANT:** The user must determine exact equilibration times for the dialysis.

4. **Change the dialysis buffer as necessary.**
5. **Pipette out the sample carefully from the *GeBAflex-tube* to a clean tube.**

If sample volume increased during dialysis, let your sample evaporate on the bench top (more details, see page 19).

### **Sample Concentration by evaporation with Mini *GeBAflex-tube***

*GeBAflex-tubes* are ideally suited for sample concentration via evaporation because of their dual membranes and large surface area. Dialysis and concentration in the same device reduce protein loss. Unlike closed-system centrifuge-type devices, sample concentration can be easily monitored in the *GeBAflex-tubes*.

1. **Place a sample in the *GeBAflex-tube* or use already dialyzed sample and place it on micro tube rack stand.**
2. **Let your sample evaporate on the bench top (using a fan to increase airflow across the membrane will speed up the process), making sure to check every 10 min or less to prevent evaporation to dryness.**

**IMPORTANT:** When evaporating water from your sample, small molecules (buffer salts, reducing agents, etc.) will also be concentrated.

### **DNA and RNA Extraction from Gel with Mini *GeBAflex-tube***

This procedure is designed to extract DNA or RNA from polyacrylamide or agarose gels.

#### **Procedure**

1. **Fill the *GeBAflex-tube* with 250 µl of dH<sub>2</sub>O, incubate for at least 5 min, empty the tube.**

**IMPORTANT:** Check carefully that no dH<sub>2</sub>O is leaking from the tube. Absorbent of water, by the dry membrane, cause to the decrease in water level.

- 2. Excise the slice of gel containing the desirable DNA or RNA fragment with a clean, sharp scalpel.**

Minimize the size of the gel slice by removing extra gel.  
Maximum gel slice size 0.4 cm x 1.1 cm.

- 3. Transfer the gel slice to a *GeBAflex-tube*. Fill the tube with 250 µl dH<sub>2</sub>O.**

Avoid air bubbles in the tube. **Don't fill the tube with several gel slices**, for larger gel slices use more than one tube.

- 4. Place the *GeBAflex-tube* in the provided tray (see Figure 1 page 8).**

The supporting tray can comprise 1-4 *GeBAflex-tube(s)*.

**IMPORTANT:** The two membranes of the *GeBAflex-tube* must be in perpendicular to the electric field to permit the electric current to pass through the tube.

- 5. Place the supporting tray containing the *GeBAflex-tube(s)* in a horizontal electrophoresis tank containing running buffer (see Figure 2 page 9).**

**IMPORTANT:** Immerse **fully** the *GeBAflex-tube(s)* with the tray in the buffer.

- 6. Pass electric current (usually at 100 volt) until the nucleic acid exits from the gel slice (see Tables 2 and 3 page 22).**

Optional: Follow the DNA or RNA eluted out of the gel with a hand-held UV lamp or table.

**IMPORTANT:** The electro-elution time need to be adjusted for each individual sample.

- 7. Reverse the polarity of the current for 120 seconds.**

This step will release the nucleic acid from the membrane.

- 8. Open the *GeBAflex-tube* gently, pipetting the solution up and down carefully (at least 5 times) and transfer the solution to a clean tube.**

Do the pipetting on the inner side of the membrane.

**Note:** Concentrate the extracted nucleic acid by standard concentration methods; for nucleic acid precipitation see page 23.

### Elution Time Tables

In this method the elution time depends on the size of the nucleic acid fragment, the concentration of the gel, the size of the gel slice, the ratio of the polyacrylamide:bisacrylamide and the applied voltage.

**IMPORTANT:** The electro-elution time at the elution step needs to be adjusted for each individual sample.

**Table 2:** Minimum time needed to extract various DNA fragments from 4% polyacrylamide gel (37.5:1 polyacrylamide:bisacrylamide) at 100 volt in 1XTBE buffer.

Fragment size (bp)	Elution time (Min)
100	10-15
200	15-20
500	30-35
1000	55-60
1400	75-80

**Table 3:** Minimum time needed to extract DNA fragments from 1% agarose gel at 100 volt in 1XTAE buffer.

Fragment size (bp)	Elution time (Min)
-----------------------	-----------------------

500	10-15
1000	15-20
2000	25-30
5000	40-45
8000	50-55
10000	55-60

## DNA or RNA precipitation

### Procedure

1. Complete the volume of the sample to 300  $\mu$ l with water.
2. Add 0.1 volume of 3M KAc pH 5.2 and 0.7-1 volume of isopropanol to the solution. Mix gently by inverting the tube several times.

For example, add 30  $\mu$ l of 3M KAc pH 5.2 and 231-330  $\mu$ l isopropanol to a 300  $\mu$ l sample.

**Note:** addition of carrier (e.g. 20  $\mu$ g tRNA or 20  $\mu$ g glycogen) to the solution will increase the efficiency of precipitation.

3. Incubate at -20°C for 10 min.

To increase DNA or RNA precipitation yield incubate the samples over night at -20°C.

4. Centrifuge the sample at 4°C for 30 min at 14,000 RPM.
5. Carefully discard the supernatant without disturbing the pellet.
6. Wash the pellet with a cooled 70% ethanol.
7. Air-dry the pellet for 5-20 min.

Do not over-dry the pellet (e.g., by using a vacuum evaporator), as this will make the DNA, especially if it is of high molecular weight, difficult to redissolve.

8. Redissolve the DNA or RNA in a suitable buffer.

Use a buffer with pH  $\geq 8.0$  for redissolving, as DNA does not dissolve readily in acidic buffers.

### Troubleshooting Guide

Problem	Cause	Comments and Suggestion
Low yield	Insufficient elution time	Increase elution time. Increase applied voltage.
	Current polarity was not reversed	Reverse the polarity of the current for 120 second.
	Incomplete emptying of the tube from the macromolecules-containing solution	Make sure to empty all the macromolecules containing solution at the end of elution.
	Ineffective precipitation	Use suitable precipitation procedures.
	Tube not fully immersed in the buffer of the electrophoresis tank	Fully immerse the tube in the buffer of the electrophoresis tank.
	Gel slice not fully immersed in the buffer inside the tube	Fully immerse the gel slice in buffer inside the tube.
	More than one gel slice was inserted into the tube	Don't fill the tube with several gel slices, for large gel slices use more than one tube.
	The electric current don't pass through the tube	The two membranes of the <i>GeBAflex-tube</i> must be parallel to the electric field.

Long elution time	Low applied voltage	Increase applied voltage.
	Gel slice is not fully immersed in the buffer inside the tube	Fully immerse the gel slice in buffer inside the tube.
	<b>Cause</b>	<b>Comments and Suggestion</b>
	Tube not fully immersed in the buffer of the electrophoresis tank	Fully immerse the tube in the buffer of the electrophoresis tank.
Macromolecules containing solution reduced after elution	Membrane not wetted before elution	Wet the membrane for 5 min with dH <sub>2</sub> O before elution.
	Pinhole in the membrane, due to careless handling of the tube	Change tube.
Presence of air bubbles in the tube	Insufficient dH <sub>2</sub> O or running buffer inside the tube	After inserting the gel slice in the tube, fill the tube to the top.

## Ordering Information

Product	Contents	Cat. No.
<b><i>GeBAflex-tube</i> kits for extraction and dialysis</b>		
Mini GeBAflex-tube (10)	10 <i>GeBAflex-tube</i> of 6000-8000 cut-off, buffers, supporting tray, floating rack	T070-6-10
Mini GeBAflex-tube (30)	30 <i>GeBAflex-tube</i> of 6000-8000 cut-off, buffers, supporting tray, floating rack	T070-6-30
Mini GeBAflex-tube (10)	10 <i>GeBAflex-tube</i> of 12000-14000 cut-off, buffers, supporting tray, floating rack	T070-12-10
Mini GeBAflex-tube (30)	10 <i>GeBAflex-tube</i> of 12000-14000 cut-off, buffers, supporting tray, floating rack	T070-12-30

Maxi <i>GeBAflex-tube</i> (5)	5 <i>GeBAflex-tube</i> of 3500 cut-off, supporting tray, floating rack	T030
Maxi <i>GeBAflex-tube</i> (15)	15 <i>GeBAflex-tube</i> of 3500 cut-off, supporting tray, floating rack	T035
Maxi <i>GeBAflex-tube</i> (5)	5 <i>GeBAflex-tube</i> of 6000-8000 cut-off, supporting tray, floating rack	T040
Maxi <i>GeBAflex-tube</i> (15)	15 <i>GeBAflex-tube</i> of 6000-8000 cut-off, supporting tray, floating rack	T045
Maxi <i>GeBAflex-tube</i> (5)	5 <i>GeBAflex-tube</i> of 12000-14000 cut-off, supporting tray, floating rack	T050
Maxi <i>GeBAflex-tube</i> (15)	15 <i>GeBAflex-tube</i> of 12000-14000 cut-off, supporting tray, floating rack	T055
<i>Midi GeBAflex-tube</i> (15)	15 <i>GeBAflex-tube</i> , 3500 cut-off, buffers, supporting tray, floating rack	T011
<i>Midi GeBAflex-tube</i> (15)	15 <i>GeBAflex-tube</i> of 6000-8000 cut-off, Buffers, supporting tray, floating rack	T021
<i>Midi GeBAflex-tube</i> (30)	30 <i>GeBAflex-tube</i> of 3500 cut-off, buffers, supporting tray, floating rack	T012
<i>Midi GeBAflex-tube</i> (30)	30 <i>GeBAflex-tube</i> of 6000-8000 cut-off, buffers, supporting tray, floating rack	T022
<b><i>GeBAflex-tube</i> kits for dialysis</b>		
Mini <i>GeBAflex-tube</i> (10)	10 <i>GeBAflex-tube</i> of 6000-8000 cut-off, floating rack	D070-6-10
Mini <i>GeBAflex-tube</i> (30)	30 <i>GeBAflex-tube</i> of 6000-8000 cut-off, floating rack	D070-6-30
Mini <i>GeBAflex-tube</i> (10)	10 <i>GeBAflex-tube</i> of 12000-14000 cut-off, floating rack	D070-12-10
Mini <i>GeBAflex-tube</i> (30)	10 <i>GeBAflex-tube</i> of 12000-14000 cut-off, floating rack	D070-12-30

Maxi <i>GeBAflex-tube</i> (5)	5 <i>GeBAflex-tube</i> of 3500 cut-off, floating rack	D030
Maxi <i>GeBAflex-tube</i> (15)	15 <i>GeBAflex-tube</i> of 3500 cut-off, floating rack	D035
Maxi <i>GeBAflex-tube</i> (5)	5 <i>GeBAflex-tube</i> of 6000-8000 cut-off, floating rack	D040
Maxi <i>GeBAflex-tube</i> (15)	15 <i>GeBAflex-tube</i> of 6000-8000 cut-off, floating rack	D045
Maxi <i>GeBAflex-tube</i> (5)	5 <i>GeBAflex-tube</i> of 12000-14000 cut-off, floating rack	D050
Maxi <i>GeBAflex-tube</i> (15)	15 <i>GeBAflex-tube</i> of 12000-14000 cut-off, floating rack	D055
Midi <i>GeBAflex-tube</i> (10)	10 <i>GeBAflex-tube</i> of 3500 cut-off, floating rack	D010
Midi <i>GeBAflex-tube</i> (30)	30 <i>GeBAflex-tube</i> of 3500 cut-off, floating rack	D012
Midi <i>GeBAflex-tube</i> (10)	5 <i>GeBAflex-tube</i> of 6000-8000 cut-off, floating rack	D020
Midi <i>GeBAflex-tube</i> (30)	15 <i>GeBAflex-tube</i> of 6000-8000 cut-off, floating rack	D022
<b>Protein Concentration kits and SDS Removing Buffer</b>		
<i>ProteoConN</i> (10)	10 <i>ProteoConN</i> columns, buffers	PN010
<i>ProteoConN</i> (20)	20 <i>ProteoConN</i> columns, buffers	PN020
<i>ProteoConN</i> (10)	10 <i>ProteoConN</i> columns, buffers	PD010
<i>ProteoConN</i> (20)	20 <i>ProteoConN</i> columns, buffers	PD020

SDS Removing Buffer	30 ml of SDS removing buffer	PDS010
SDS Removing Buffer	60 ml of SDS removing buffer	PDS020
<b>Protein Staining Solution</b>		
<i>SeeBand</i> protein staining solution	500 ml solution and handbook	SB010
<i>SeeBand Forte</i> protein staining solution	500 ml solution and handbook	SB020